



902 13th Street North  
 Benson, MN 56215  
 320-843-4109  
 Fax 320-843-2074

Highway 15  
 Northwood, ND 56267  
 701-587-6010  
 Fax 701-587-6013

## MANURE INFORMATION FORM

Fill bottle 3/4 full. DO NOT OVERFILL. Close cap securely. Write sample ID on bottle cap. Put bottle in plastic bag and seal.  
 Completely fill out information Form. PLEASE DO NOT PUT FORM IN PLASTIC BAG WITH SAMPLE.  
 To preserve sample during shipping, freeze or use ice pack.

### Grower:

Name
Address
Address
City, State, Zip

### Submitter:

Account Number
Name
Address
City, State, Zip

### SAMPLE IDENTIFICATION:

<p style="text-align: center;">Test Analysis</p> <p><b>Complete:</b> % Moisture, Total Nitrogen, Phosphorus, Potassium, Sulfur, Calcium, Magnesium, Sodium Manganese, Copper, Zinc, and Iron <input type="checkbox"/></p> <p><b>Complete with Ammonia:</b> <input type="checkbox"/></p> <p><b>Partial:</b> % Moisture, Total Nitrogen, Phosphorus, Potassium <input type="checkbox"/></p> <p><b>Partial with Inorganic N:</b> % Moisture, Total Nitrogen Phosphorus, Potassium, Ammonia, Nitrate <input type="checkbox"/></p> <p><b>Partial with Ammonia:</b> <input type="checkbox"/></p>	<p><b>ADDITIONAL TESTS:</b></p> <p>Ammonical Nitrogen <input type="checkbox"/></p> <p>Nitrate Nitrogen <input type="checkbox"/></p> <p>Total Carbon <input type="checkbox"/></p> <p>pH <input type="checkbox"/></p> <p>Salts <input type="checkbox"/></p> <p>Volatile Solids <input type="checkbox"/></p> <p>Percent Ash <input type="checkbox"/></p>	<p><b>ADDITIONAL TESTS:</b></p> <p>Sulfur <input type="checkbox"/></p> <p>Zinc <input type="checkbox"/></p> <p>Iron <input type="checkbox"/></p> <p>Copper <input type="checkbox"/></p> <p>Manganese <input type="checkbox"/></p> <p>Sodium <input type="checkbox"/></p> <p>Calcium <input type="checkbox"/></p> <p>Magnesium <input type="checkbox"/></p> <p>Chloride <input type="checkbox"/></p>
--	---	---

Sample Source		Type		Storage	
Dairy <input type="checkbox"/>	Poultry <input type="checkbox"/>	Solid <input type="checkbox"/>	Lagoon <input type="checkbox"/>	Slurry <input type="checkbox"/>	Pit <input type="checkbox"/>
Calf <input type="checkbox"/>	Beef <input type="checkbox"/>				
Swine <input type="checkbox"/>	Other <input type="checkbox"/>	Liquid <input type="checkbox"/>	Other <input type="checkbox"/>		

LAB USE ONLY

### SAMPLE IDENTIFICATION:

<p style="text-align: center;">Test Analysis</p> <p><b>Complete:</b> % Moisture, Total Nitrogen, Phosphorus, Potassium, Sulfur, Calcium, Magnesium, Sodium Manganese, Copper, Zinc, and Iron <input type="checkbox"/></p> <p><b>Complete with Ammonia:</b> <input type="checkbox"/></p> <p><b>Partial:</b> % Moisture, Total Nitrogen, Phosphorus, Potassium <input type="checkbox"/></p> <p><b>Partial with Inorganic N:</b> % Moisture, Total Nitrogen Phosphorus, Potassium, Ammonia, Nitrate <input type="checkbox"/></p> <p><b>Partial with Ammonia:</b> <input type="checkbox"/></p>	<p><b>ADDITIONAL TESTS:</b></p> <p>Ammonical Nitrogen <input type="checkbox"/></p> <p>Nitrate Nitrogen <input type="checkbox"/></p> <p>Total Carbon <input type="checkbox"/></p> <p>pH <input type="checkbox"/></p> <p>Salts <input type="checkbox"/></p> <p>Volatile Solids <input type="checkbox"/></p> <p>Percent Ash <input type="checkbox"/></p>	<p><b>ADDITIONAL TESTS:</b></p> <p>Sulfur <input type="checkbox"/></p> <p>Zinc <input type="checkbox"/></p> <p>Iron <input type="checkbox"/></p> <p>Copper <input type="checkbox"/></p> <p>Manganese <input type="checkbox"/></p> <p>Sodium <input type="checkbox"/></p> <p>Calcium <input type="checkbox"/></p> <p>Magnesium <input type="checkbox"/></p> <p>Chloride <input type="checkbox"/></p>
--	---	---

Sample Source		Type		Storage	
Dairy <input type="checkbox"/>	Poultry <input type="checkbox"/>	Solid <input type="checkbox"/>	Lagoon <input type="checkbox"/>	Slurry <input type="checkbox"/>	Pit <input type="checkbox"/>
Calf <input type="checkbox"/>	Beef <input type="checkbox"/>				
Swine <input type="checkbox"/>	Other <input type="checkbox"/>	Liquid <input type="checkbox"/>	Other <input type="checkbox"/>		

LAB USE ONLY